

On September 28, 2020, ECG Management will Go-Live to increase the efficiency in which ECG tracings are imported into Cerner. The preferred workflow is to start with an order however there are times when the patient is not registered, and an ECG must be obtained emergently. This flyer outlines the steps for **Clinical/Tech** staff to match an ECG tracing to an ECG or Stress Test order enabling the study to be electronically available for the Cardiologist to read in ECG Management.


Unmatched Study

ECG orders should be entered before the exam is acquired for it to be electronically reviewed by a Provider. ECGs should only be printed on an as-needed basis, such as when a patient has **Not** been registered in Cerner and needs an immediate ECG printout. If an ECG is acquired prior to order entry the study imports into ECG Management in an **Unmatched** state. The study must be matched to the order for it to be electronically available in the patient chart for the ordering provider to review and for the Cardiologist to perform their final interpretation.

Unmatched Study



If an ECG study is acquired without an order the staff member who obtained the ECG is responsible for entering the appropriate order in the patient chart and matching the order to the ECG tracing in ECG Management.

STEP 1: From the patient chart enter the appropriate ECG or stress test order.

STEP 2: Click the Cardiovascular button in the toolbar. 


STEP 3: Locate the appropriate patient within the Procedures worklist using the **column filter**.

- Right click patient row
- Hover over **Action**
- Hover over **Match Study**
- Click **Match Study**
- **No Results Found** window will pop up, click **OK**.




STEP 4: Match Study Window

- Uncheck **Filter by Last Name**
- Select **Modify Search**




STEP 5: Unmatched Studies window

- Use the **Date Range** dropdown to select a time range.
- Click the check box to the left of ECG
- Click **OK**



STEP 6: Match Studies window

- Locate the appropriate patient by comparing the top **Procedure to Match** pane to the bottom **Select an Unmatched Study** pane.
- Click the appropriate patient row to highlight
- Click **Match**



STEP 7: Study State updates to Matched Viewable.

- The study is now available to the Cardiologist to complete their final interpretation.

Cardiovascular				
Last Name	First Name	Type	Status	Study State
TESTING	CIECG73	CV ECG	Ordered	Matched Viewable

NOTE: ECGs tracings cannot be matched to orders on the ECG carts. Matching can only occur within ECG Management.


Unmatched Study without Order

At the end of every shift each department should check the **Unmatched Study List** and the ECG carts to ensure all tests have orders and were transferred into Cerner. ECG tracings that populate the Unmatched Studies list may not have ECG or Stress orders and will need to have the appropriate order entered before the tracing can be imported to Cerner.

STEP 1: Click the **Cardiovascular** button in the toolbar.




STEP 2: Click **Cardiovascular** in the top toolbar



- Select Show **Unmatched Studies...**


STEP 3: **Unmatched Studies** window

- Use the **Date Range** dropdown to select time range.
- Click check box to the left of ECG
- Click **OK**



STEP 4: **Unmatched Studies** list

- The unmatched studies list defaults to all member organizations.
- Click the **My Orgs** dropdown to filter by your location.



NOTE: Location will not save as a preference for the Unmatched Studies list and will need to be selected every time the list is reviewed.

STEP 5: Navigate to the patient's chart using the **FIN#** associated with the ECG acquisition.

STEP 6: Enter the appropriate order.

STEP 7: Click the **Cardiovascular** button to navigate back to ECG Management.

STEP 8: The study updates on the Procedures worklist to an **Ordered, Unmatched** state.

STEP 9: Follow the **Unmatched Study** steps above for matching orders to ECG tracings.

NOTE: Click [here](#) to watch a Cerner developed education simulation on Matching an ECG to an Order.