


Agfa EI is an application where cardiac images are stored, displayed, and reviewed.

Agfa EI Icon Tips and Tricks



Cerner/Epic and Agfa EI are integrated.





- When a patient is opened in Cerner or Epic, Agfa EI opens to the image area.
- Images are displayed, and enhancement tools can be used for review.







Change the layout to increase or decrease how many images appear in the viewport area by clicking the 1 to change the selection.

Hot Spot Tools






➤ Left Side of the Image Area






	Zoom Tool
	Pan Tool
	Window Leveling
	Flips image horizontally, vertically, and rotates 90 degrees clockwise.

➤ **Right Side of Image Area**

	<ul style="list-style-type: none"> • If the patient has prior studies in the database, they are shown in the patient history tab. • Click a Prior Study in Patient History to see the thumbnail previews. <ul style="list-style-type: none"> ○ Images of a Prior Study are surrounded with light blue frames, including study date and time. ○ Drag and drop one of the thumbnails into the workspace in order to compare images of the primary study with images of the prior study.
	<p>Review icon – this is where shutter or window leveling is used.</p>
	<p>Measurement icon – used to make measurements.</p>
	<p>Tool icon – use to do annotations and anonymize.</p>

➤ **Additional Icon Actions**

	<p>Left click into the image or cine loop to insert free text on one selected frame.</p>
	<p>Insert an arrow at one selected frame of a cine loop.</p>
	<p>Removes the last modified annotation of the selected cine loop or image.</p>
	<p>Increase or decrease playback speed.</p> <ul style="list-style-type: none"> • Range is 10-100% and will change in 10% increments.
	<p>Delete all Modifications – Remove all annotations of the selected cine loop or image.</p>

	<p>Export the complete content of the review area by clicking the Save AVI or BMP button.</p> <ul style="list-style-type: none"> • If at least one loop is playing, the complete screen is exported as AVI, otherwise as BMP, PNG, or JPG.
	<p>Create a PDF document including all the images which are shown in the workspace.</p> <ul style="list-style-type: none"> • For more images, change the tiling of the workspace.
	<p>Bottom of Image Area</p> <ul style="list-style-type: none"> • Use these tools to move through images in series and cines.
	<p>Reset – Displays the dataset as originally displayed in the beginning of each workflow step.</p>
	<p>Save Secondary Capture – Store results on the images and cine loops using this button.</p> <ul style="list-style-type: none"> • The complete image region is captured. • The new object becomes visible in the thumbnail preview.