
XERO Viewer is a medical image viewer that can be launched by an image link in the medical record or from a Web-Browser

Accessing XERO Viewer

- Open the patient chart in PowerChart, navigate to the **Imaging component** on a **Workflow MPage**, and select the imaging study to be viewed.

- Click the **View Image** icon. 

NOTE:  **Image link** in the EHR (Electronic Health Record) opens directly to the patient and study that is being opened. The image status must be Final Report to view an image from the EHR.

Accessing XERO Viewer using a Browser

- Users not needing to use Ascend reporting.
 - Click to open link <https://brnagfaweb.me.emh.org/>
- Cardiology and EMMC Vascular workflow to use Ascend reporting.
 - Click to open link <https://brnagfaweb.me.emh.org/?returnTo=/workflow/>

NOTE: These links can be added to your **My Quick Links** on the **Intranet Home Page**.

STEP 1: Enter **Patient Name**, **Date Of Birth**, **Accession number**, and **Study Date** (if known), and **Modality** to search for patient.


STEP 2: Click **Search**.

STEP 3: Select method of display, **Report and Images**, **Report Only**, or **Image Only**.



STEP 4: Double-click the study.

STEP 5: Select right arrow to show image area.



From the Office of Clinical Informatics

Xero Viewer

March 27, 2024

Page 2 of 2

STEP 6: Select arrow on the bottom middle of the screen to scroll between images.



STEP 7: **Patient Banner** is used to verify patient.



STEP 8: Select the switch series icon in the upper right to switch between series.

STEP 9: Click one of the **Tools** to assist with reviewing images.



STEP 10: Select the ellipsis to change layout display.



STEP 11: Use the **Help** link to assist with any questions.

STEP 12: Click **Logout** once review of image(s) has been completed.

