
The Cerner Patient Observer (CPO) solution allows NLH to virtually monitor patients identified as a high-fall risk using 3D cameras and sophisticated algorithms to track patient movement. Motion zones are defined to meet the unique needs of each patient and two-way communication helps virtual observation staff redirect the patient and potentially prevent an adverse event.

Onboarding

- Camera is seen but offline.
 - Instruct the nurse to unplug the camera from the power pack and plug back in.
 - Nurse should hear several beeps from camera when it is plugged back in.
- Cannot See Patient.
 - Ask nurse to remain in the room.
 - Click **Reboot Camera**.
 - If not successful, click **Secondary tab** and select **Reboot Room PC**.
 - Ask for a different camera to be brought to the patient room if none of the above work and ask nursing to place a HelpDesk ticket to have the camera fixed.

Motion Zone


- Motion zones too sensitive and are triggering false movement.
 - Click **Edit Patient Zones**.
 - Click **plus sign icon**. 
 - Double-click inside rectangle.
 - Increase **Min Blob Size** by increments of 25 to decrease the motion sensitivity.
 - Click **Save**.
 - Re-evaluate the location of the motion zones to determine if they are too close to the patient.

Mouse

- The mouse pointer cannot be located.
 - Touch the screen on the small monitor to bring the mouse pointer to that location.


Headset

- Headset is not working, not able to hear.
 - Click the **Microphone** button in the upper left side of the large monitor.
 - If it defaults to **None**, click the available **Headset** option.
 - Click **Save**.




IBUS

- The iBus (communication system that transfers data between components or computers) square is red with an X instead of being green and red eyeglasses display.
 - Call to notify the nurse video feed has been lost and someone needs to be with the patient until the issue is resolved.
 - Call the HelpDesk.



How can I tell when iView documentation was done last?


- The Last Done Date and Time column in the Multi-Patient Task List (MPTL) displays when documentation was last documented in iView.



How can I complete documentation on a patient that fell off the MPTL?

STEP 1: Select **Options** in the toolbar at the top of the screen.

- Then click **Task Display**.



STEP 2: Under **Status**, check the box next to **Discontinued/Canceled**.

STEP 3: Click **OK**.

STEP 4: Patients with orders that have been discontinued will have **red** in the first column and **Discontinued** under **Task Status**.

Location/Room/Bed	Name	Task Status	Scheduled Date and Time	Task Description	Order Details	Last Done Date and Time
SBLA / B584 / 01	*TESTING, WILLIAM 01/06/1960 <i>*Name Alert</i>	Discontin...	11/21/2022 13:15 EST	Virtual Patient Observation	Start: 11/21/22 13:15:00 EST, History of falls	
SBLA / B583 / 01	*TESTING, OLLIE 06/07/2001 <i>*Name Alert</i>	Pending	PRN	Virtual Patient Observation	Start: 11/21/22 14:22:00 EST, History of falls	
SBLA / BHAL / 03	*TESTING, ROMAN 07/23/1985 <i>*Name Alert</i>	Pending	PRN	Virtual Patient Observation	Start: 10/31/22 12:38:00 EDT, History of falls	11/21/2022 13:39 EST

STEP 5: To sort discontinued orders from active orders, double-click **Scheduled Date and Time** column.

STEP 6: To document in iView, double-click anywhere in the patient row.

STEP 7: Once the Virtual Observation Technician is done documentation, remove patients with discontinued orders by following Step 1.

- Under **Status**, uncheck the box next to **Discontinued/Canceled**.
- Click **OK**.