


This flyer will introduce you to the Patient Demographic Bar as well as the various action buttons which will assist you in navigating the Blood Transfusion Workflow.

Reviewing the Main Menu

After you identify the patient, the patient demographic information is displayed at the top of the main menu page.

➤ Patient Demographic Bar




- **ID Number:** This is the patient's Financial Identifier Number (FIN) number.
- **Patient Age:** If the patient has been marked as deceased, the age is followed by (D).
- **Patient Blood Type:** The box displays the patient's blood type that was entered during the first Start Transfusion workflow. If no blood product has been started, the value for the Blood Type box is **Unk**. After the patient has had one unit transfused and has been type and crossed and has a component tag associated, it will display the blood type. This will update when the transfusion has started.

NOTE: All Patient Information in the demographics bar, except the blood type, is displayed through an interface from the hospital's admission, discharge, and transfer (ADT) system.

➤ Transfusions

- There are multiple action buttons to assist with transfusion documentation.
 - **Browse History:** Use the Browse Transfusion History Page to review a summary of the administration activity for your patient.
 - **View Reaction:** This displays the clinical signs or symptoms that were selected from the list of hospital defined transfusion reactions during, End Transfusion, or Post Transfusion Reaction processes.




➤ Vital Signs

- **Record Vital Signs:** You will record vital signs during the Start, End and Post-Transfusion Reaction workflows. You access Record Vital Signs from the patient's Main Menu Page.
- **Browse Vital Signs History:** Use Vital Signs History Page to review all the vital signs activity recorded for your patient.



NOTE: After you complete an action, the system always displays the patient main menu. From the patient main menu, you can view the Identify Patient page or Log Out.

- To view the Identify Patient page, click the close button (X) in the upper-right corner.



- To log out, click Log Out.

