

The Clinical Leader Organizer (CLO) has icons that fire based on orders and nursing documentation within the EHR. These icons will assist leaders and nursing staff in identifying if documentation for certain measures has been completed as well as pinpointing certain risk factors.

What do the Icons indicate and how do they populate on the Clinical Leader Organizer?

- Click on an icon for more documentation or order details.
- Periodically refresh the CLO to get the most recent data. It does not automatically refresh.

NOTE: These icons may not match the icons for the same thing in Capacity Management.

ICON OR INDICATOR	DESCRIPTION	AVAILABLE COLUMN
	This icon displays when the patient has an artificial airway documented. Click the icon to see the type of artificial airway, Endotracheal Tube Activity, date and time of last assessment and who documented the assessment.	Airway
	This icon displays when the patient has documented allergies. Click the icon to view the allergen, reaction, and severity documentation.	Allergy
	This icon displays when No Allergy documentation has been recorded.	Allergy
	This icon displays when a free text allergy has been documented.	Allergy
	This icon displays when No Known Medication Allergies are documented.	Allergy
--	A dash indicates information is available. Click in the Care Team column to see who has established a relationship with the patient. Click Display Physician to have the provider name populate in the column for all patients on the list. Display Nursing becomes available to populate the column once Display Physician is selected. Note that nursing names do appear under Physicians and will only appear under Assigned Nurses if the nurse is assigned using Staff Assignment in Clairvia. Only one name will populate the column, a physician or a nurse.	Care Team

	This icon displays when the Urinary Catheter dynamic group has been created and Insert New Catheter or Daily Necessity Assessed has been documented in Activity. Click the icon to see which Activity was documented.	Catheter
	This icon displays when no documentation in Activity has been documented in the last 72 hours.	Catheter
	This icon displays when an active order for a Central Line is present. It also displays when Activity and dressing change is documented.	Central Line
	This icon displays if Central Line dressing documentation has exceeded 7 days (168 hours).	Central Line
	This icon indicates the patient has a diet order. Click the icon to view the diet and order details.	Diet
	This icon displays when the NPO Communication order is placed.	Diet
	This icon displays if the patient has more than one diet order. This includes a diet order and an NPO communication order.	Diet
	This icon displays when the patient has an active discharge order.	Discharge
	This icon displays when the patient has been discharged. Click the icon to display the date and time of discharge.	Discharge
	This icon will display if the Elopement Precautions order has been entered on an Acadia Hospital patient encounter.	Elopement
	This indicator displays the most recently documented fall risk score. A trending up or down arrow will indicate the score trend. If no change from the previous value, an arrow will not display. Click the icon to see the last assessment date and who completed it.	Falls Risk
	This icon displays when the patient has one or more of the following orders placed: Precautions (Nsg) - any of the options in the order, Continuous Observation , or Room Search (Nsg) . Click the icon to see the type of precaution & the order date.	High Risk
No Icon	The location column displays the nursing unit name, room number, and bed placement.	Location
	This icon displays when the Medication History hasn't been started.	Medication History

	This icon displays when the Medication History has been completed. Click the icon to see when it was completed and by whom.	Medication History
OBS	This indicator displays when the patients encounter status is Observation. No data displays if icon is clicked.	Observation
O2	This indicator will display if Initiate is documented in Oxygen Activity. Click this icon to see the FIO2, and Flow rate, when and by whom this information was documented.	Oxygen Therapy
7	This indicator displays most recently documented pain score. Click the icon to see the last three pain scores documented. Click Expand All to see when the score was documented and by whom.	Pain
No Icon	The Patient column displays the patient name, age and gender. Click the patient name to open the patient chart.	Patient
	This icon displays when an order for Behavioral or Non-Behavioral Restraints has been placed. Click the icon to see the order details.	Restraints
	This icon displays when the patient has a code status of Full Resuscitation . Click the icon for order details.	Resuscitation
	This icon displays when the patient has a code status of DNR (aggressive care unless restricted) , or DNR with Comfort Care Only . Click the icon for order details.	Resuscitation
	This icon displays if the patient has a code status of Full Code, but Do Not Intubate . Click the icon for order details.	Resuscitation
10 ↘	This indicator displays the most recently documented Braden Score . A trending up or down arrow displays based on a comparison of the most recent documentation. No arrow indicates the same value. Click the icon to view the last 3 documented scores. Click Expand All to view when documented and by whom. There is an Idea on the Cerner Idea Page to have the number display in red for at risk scores and is listed for Future Roadmap Consideration. Consider giving it support. https://connect.cerner.com/ideas/26059#comment-573386	Skin Integrity
	This icon displays when the patient has a Suicide Precautions order and/or if Suicide Ideation is documented. Click the icon to see the date and time ordered, and the frequency of the Suicidal Ideation, along with last assessment date & time and who documented the data.	Suicide

	This icon displays if Suicidal Ideations are documented and the patient does not have a Suicide Precautions order.	Suicide
	This icon indicates the patient has a Telemetry Monitoring order . Click to view the date and time ordered.	Telemetry
	This icon displays when there is a Vent Settings order . Click the icon to view the order details.	Ventilator
No Icon	Displays the Length of stay (LOS). Click to view the admission date. The LOS does not update if the patient changes from Observation to Inpatient status.	Visit