

#### From the Office of Clinical Informatics

# Cerner Millennium Activating Future Orders

August 11, 2023

Patients who show up at a Quest Diagnostics draw station who initially were intended to go to the Hospital or Outside Lab draw station will need to have the order activated.

Recent •

### **Searching for the Patient**

**STEP 1:** Sign into **PowerChart**.

**STEP 2**: Click the **magnifying glass**.

**STEP 3:** Using **Search**, search for the patient's name.

STEP 4: Locate the most recent practice encounter and select that FIN.

 If no encounter found, inform Quest Lab Staff to have patient contact their providers office to confirm the ordering providers location.

**STEP 5:** Click **OK** to open the patient chart.

# **Locating Ordering Provider Name**

**STEP 1:** Locate **Orders** in the blue menu to the left and right-click.

**STEP 2:** Select **Set as Default View**.

• Next time a patient's chart is entered, it will default open to Orders.

**STEP 3:** Click **Orders** link in the blue menu.

Orders Profile page opens.


**STEP 4:** Locate the intended Order from the Laboratory View.


**STEP 5:** Locate the order, right-click, and select **Order Information**.

- Locate ordering Provider.
- Close Order Information box.









**NOTE:** 

If unable to locate the order, inform the Quest Lab Staff to have patient contact their providers office to confirm the order(s) have been placed.

## **Modifying and Activating the Order**


**STEP 1:** Locate order in the Order Profile page, right-click order, and select **Modify.** 

STEP 2: Enter the **ordering provider's name** and select **Administrative** for the Communication Type.

**STEP 3:** Change **Performing Location (ARLN)** to **Quest Diagnostics**.

**STEP 4**: Click **Orders for Signature**.





STEP 5: Click Sign.

STEP 6: Right-click the order and select **Activate**.

**STEP 7:** Click **Orders for Signature**.

STEP 8: Click Sign.

