


The Quick Connect workflow allows providers to connect remotely to the Amwell carts onsite at Northern Light facilities.

### Connecting to Amwell Cart via Quick Connect

- STEP 1:** Click **Clinical Links** in the weblinks toolbar in Cerner.
- STEP 2:** Click **Amwell Telehealth Quick Connect**.
- STEP 3:** Sign in with appropriate email address and password.
- STEP 4:** Type the patient's room number in the **Room name** field.
- STEP 5:** Click **Select invitation method** and select **Endpoints**.
- To connect to a **laptop via an email**, select **Email** for the invitation method and enter the appropriate email and navigate to Step 8.
- STEP 6:** Click **All Facilities** and select the appropriate facility.
- STEP 7:** Click **Search** and select the appropriate cart.
- NOTE:** **Communication to the connecting site may be required to ensure the correct cart is called.**
- STEP 8:** Click **Invite Participants**.
- This initiates the call into the selected cart.
- STEP 9:** Click **Join Room Now** to join the call.
- STEP 10:** At the end of the call, refer to **Recent Connections** as needed.
- NOTE:** **If needed to reconnect to a cart, send a new invite to the cart by completing steps 5-9.**



### Downtime Procedure

During a Cerner downtime, this workflow remains the same.